ERA LUMINA

BiOTEK
MILANO

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V.01.2024

1. Introduction	03
2. General warnings	03
3. General notes on delivery	04
4. Fields of application	04
5. Safety	05
6. Technical safety	05
7. Safety of use	05
8. Accessory safety	06
9. Information for waste disposal	06
10. ERA handpiece diagram	06
11. ERA handpiece technical data	07
12. Operational connections and assembly	07
13. Electrical connection	08
14. Preparation of work area and device	08
15. Installing the cartridge	08
16. Needle adjustment	09
17. Cartridges available	09
18. General advice of use	09
19. Troubleshooting	09
20. Maintenance and cleaning	10
21. Declaration of conformity	11
22. Warranty certificate	11
23. Returns	12
24. Disputes	12

Dear Customer,

Thank you for purchasing this BIOTEK product, we are sure that you will be satisfied with it.

ERA is an original product, designed and manufactured at the BIOTEK facility in Italy.

BIOTEK has been producing device and accessories for permanent make-up, professional dermopigmentation and tattooing since 1985.

Many years of experience ensures the quality of a product specifically designed for the unique needs of the professional.

By selecting the best materials and components, BIOTEK guarantees safety, reliability, precision and high performance.

The device, subject of this manual, complies with the regulations in force in Europe and the United States. If used correctly according to the manufacturer's instructions, it does not present any danger for the operator or the customer.

Consult the website www.biotekmilano.com for further information and specifications.

The instructions for use are valid for ERA device and related accessories and contain important information for commissioning, operating and maintaining the product safely and in accordance with the intended use

IMPORTANT

THIS MANUAL REPRESENTS AN ESSENTIAL ELEMENT FOR THE CORRECT AND SAFE USE OF ERA DEVICE AND RELATED ACCESSORIES

- Read the following carefully and observe the instructions before operating the device.
- Keep these instructions for use so that they are accessible at any time to all people responsible for using, cleaning and storing the device.
- Pay particular attention to the chapters marked with the symbol: 

1. INTRODUCTION

BIOTEK device was designed by our engineers and technicians in collaboration with the major specialists in the field.

All devices are subjected to verification tests at each production stage and tested individually before being introduced onto the market in accordance with the internal production protocol.

BIOTEK device meets professional needs and is therefore suitable for those operators with specific operational skills, as well as authorized for use by the laws in force in different countries.

2. GENERAL WARNINGS

This manual is intended to provide instructions related to:

- correct assembly of components;
- safe operation;
- regular maintenance;
- information related to risks for the patient and user.

The device must be used in compliance with the procedures contained in this manual and never for purposes other than those foreseen therein..

ATTENTION

The manufacturer cannot be held responsible for incorrect use of the device or for its unsuitability for use by the purchaser who must inform himself about the laws in force in his country for carrying out the treatment and obtain himself with any professional certificates or licenses to practice. BIOTEK declines any responsibility for incorrect use of the product, for lack of knowledge of the methods of use, as well as for the non-application of procedures and regulations in the hygiene and health sector.

 ***The product must be connected to a power supply network designed in accordance with current regulations. Use, cleaning and maintenance are prohibited for minors unless assisted by an expert adult who assumes responsibility. Never leave the device or its parts within the reach of children, they could cause serious injuries. The device is intended for professional use and as such must only be used by those who hold this qualification.***

3. GENERAL NOTES ON DELIVERY

Upon receipt of the product, check that the packaging is intact.
Open the package and check that:

- the supply corresponds to the specifications of your order;
- the product has no apparent damage;
- the product contains all pieces relating to the type ordered, according to the list shown on the BIOTEK price list in force on the date of the order.

In case of damage or missing parts, immediately inform your dealer and/or BIOTEK srl.

! IMPORTANT

Carefully check that your product contains the references to access the usage and maintenance manual, the CE certificate and the warranty.

4. FIELDS OF APPLICATION

The **ERA** device, with its disposable accessories and specific pigments, allows you to successfully practice both permanent make-up, semi-permanent make-up, aesthetic and reconstructive dermopigmentation, as well as artistic tattoos of small dimensions through the introduction of pigments under the skin conveyed by a non-hollow needle inserted into a disposable hygienic module (cartridge).

Main fields of application:

- **aesthetic micropigmentation, semi-permanent make-up, permanent make-up**
- **corrective aesthetic tattoo or dermal design**
- **dermography and reconstructive dermopigmentation**
- **scalp pigmentation**
- **tattoo**

! *The safety of the product is guaranteed not only by the use of this manual, but also by the following rules:*

- ***The product must be connected to a power supply network whose system has been built according to current system standards.***
- ***Use is prohibited for minors unless assisted by an expert adult who assumes responsibility.***
- ***Never leave the device or its parts within the reach of children, as it can cause serious injuries.***
- ***The device is intended for professional use and as such must only be used by those with suitable qualifications.***

5. SAFETY

-  **ERA** device is not designed for use in the presence of explosive gases or vapors.
-  Do not pour water or liquids of various kinds on the device to avoid causing a short circuit and damaging the external coatings with corrosion processes.
-  Before any maintenance or cleaning of the device it is absolutely necessary to remove the electrical voltage from the power supply by unplugging the plug. **DO NOT** use autoclaves or ultrasonic cleaners to clean the handpieces, console and related accessories. **NEVER** immerse in water or disinfectant solution.
-  Never open the device for any reason.
-  Non aprire mai per nessun motivo l'apparecchiatura.
-  Request only original BIOTEK spare parts.

ATTENTION

The device can only be used for carrying out procedures foreseen in the specified fields of application. Other uses are at your own risk and may not be tolerated by the device. The manufacturer is not responsible for any damage caused by improper use of the device.

Use only on dermal surfaces with the exception of orifices, mucous membranes, eyes and damaged skin or near moles or skin irritations. The device cannot be used in the presence of relative humidity higher than 60%.

ATTENTION

Do not allow children or people who do not have specific skills to use the device or its accessories for any reason.

6. TECHNICAL SAFETY

- The device can be used with BIOTEK Console (mod. B4U) or with a different power supply with voltage 5.0 – 8.0 VOLT DC.
- Before using the **ERA** device, check that the voltage values shown on the power supply are suitable for the handpiece and that the polarity is correct (central pole connected to +V).
- Do not use the device if the power cord is damaged or defective.
- Do not use the device if you notice that the cartridge assembly is not locked.
- A damaged component must be replaced as soon as possible at the BIOTEK assistance center in Milan.
- Any repairs can only be carried out by BIOTEK technicians or by personnel authorized by BIOTEK.
- Repairs that are not authorized or carried out by non-BIOTEK personnel may void the safety of the device.
- Unauthorized repairs will invalidate the warranty.

7. SAFETY OF USE

- The safe use of the device cannot leave aside in-depth knowledge of adequate operating techniques and the implementation of all hygienic-sanitary precautions.
- Use only original BIOTEK or compatible sterile cartridges.
- The cartridges are marketed in sterile packaging, check the sterility expiry date shown on the label. **NEVER** use expired cartridges. Check the integrity of the packaging before use.
- **NEVER** start the device with the needle resting on a body surface.

8. ACCESSORY SAFETY

Cartridges are disposable accessories, they cannot be reused.

The cartridge hygiene module (single needle/tube system) is disposable. Do not use the same one for more than one procedure and for multiple clients. Each treated customer is entitled to receive a new set of disposable accessories. Always maintain maximum hygiene and safety in the use of your device.

Use original BIOTEK accessories or compatible cartridges.

9. INFORMATION FOR WASTE DISPOSAL

The cartridges act as a needle disposal container and can be disposed of without further protection, however they are considered special waste and must be disposed properly. Find out about the laws in force in different countries. In any case, never throw needles or cartridges in normal urban or domestic waste.

10. ERA LUMINA HANDPIECE DIAGRAM

Parts descriptions

1. Color tank
2. Disposable needle cartridge
3. Ring for needle exit adjustment
4. Handpiece body
5. Cable Connector

11. ERA LUMINA HANDPIECE TECHNICAL DATA

TYPE OF DEVICE: HANDPIECE FOR TREATMENTS
MANUFACTURER BIOTEK SRL VIA FARNETI 8 20129 MILAN ITALY
POWER SUPPLY 5.0 - 8.0 VOLT DC
FUSE NOT NECESSARY
ABSORBED POWER 1.6 W UNDER LOAD
DIAMETER 23.5 mm.
LENGTH 113 mm.
WEIGHT 110 gr.
TYPE OF CE STANDARDS
HIGH PERFORMANCE CORELESS MOTOR
N° OF NEEDLE STROKES: 85 TO 140 PER SECOND

12. OPERATIONAL CONNECTIONS AND ASSEMBLY

The **ERA** device is tried and tested directly in the factory and does not require any further calibration before installation. For commissioning, connect to a compatible power supply or to a BIOTEK Console (mod. B4U). Take the supplied connection cable and insert the connector (A) into the handpiece outlet (B). The connector is inserted with a simple but firm pressure. Insert the other end of the cable (C) into the output of the BIOTEK Console (mod. B4U) or an alternative power supply compatible with an RCA connector.

Now the **ERA** system is ready to be turned on and adjusted, depending on the type of work to be performed.

13. ELECTRICAL CONNECTION

Before connecting the power supply to the mains, take the following precautions:

1. Place the **ERA** device on a flat, stable surface.
1. Do not install the power supply and **ERA** device in excessively humid environments.


14. PREPARATION OF WORK AREA AND EQUIPMENT

Once all the required connections have been made, before starting the treatment, set up the work area as required by the hygiene protocol. The assembly of sterile and disposable accessories (cartridge) and their adjustment must only take place close to the treatment itself.

15. INSTALLING THE CARTRIDGE

Parts descriptions

1. Needle
2. Cartridge
3. **ERA** handpiece
4. Needle exit adjustment

To carry out a procedure it is necessary to use an original BIOTEK or compatible accessory (cartridge). The cartridge consists of a single unit consisting of a tip or tube which contains the single needle or several assembled needles inside. A special protection membrane is inserted into the cartridge, mechanically preventing the passage of color from the cartridge to the handpiece. The whole unit constitutes a disposable set supplied in sterile packaging. The cartridge system guarantees ease of assembly, working precision and meets the standards defined in the operational hygiene protocol.

Follow the assembly and use instructions to avoid malfunctions or failures. The cartridge (2) must be inserted on the terminal of the ERA handpiece (3). The assembly operation is very simple to perform.

Assemble and prepare as follows:

Put on disposable gloves and remove the cartridge by opening the sterile package. Insert the cartridge on the handpiece into the appropriate hole (3). The cartridge should be rotated clockwise a quarter turn until it locks and remain secure on the handpiece.

The system is now ready for use: by activating the device from the console, the needle (1) included into the cartridge (2) will automatically come out. Proceed with the adjustment of the needle outlet (4) if necessary (chap. 16). At the end of the work, with the machine turned off, it is possible to remove the cartridge from the handpiece by turning it anti-clockwise and throw it in special waste. Always use gloves and pay attention to accidental puncture.

**! IMPORTANT: NEVER GRASP THE CARTRIDGE BY THE TIP
ALWAYS TURN OFF THE DEVICE BEFORE FITTING OR REPLACING A CARTRIDGE**

16. NEEDLE ADJUSTMENT

ERA is equipped with micro-adjustment of the needle exit.

The needle exit from the cartridge is already factory calibrated at 2.5 mm

If a micro-adjustment is necessary, this can be done by rotating the appropriate ring (chap. 15, point 4).

Attention. ADJUSTMENT OF THE NEEDLE EXIT MUST ALWAYS BE CARRIED OUT WITH THE DEVICE TURNED ON. WITH THE DEVICE TURNED OFF IT CANNOT BE DETERMINED EXACTLY HOW MUCH THE NEEDLE WILL PROTRUDE ONCE IT IS ACTIVATED. NEVER ROTATE THE CARTRIDGE TO ADJUST THE NEEDLE OUTPUT AS THE CARTRIDGE COULD BE DISCONNECTED FROM THE HANDPIECE

17. CARTRIDGES AVAILABLE

The operator will be able to choose from a wide variety of types of tips available depending on the type of treatment to be performed and the operating method. Consult the current price list to view the BIOTEK range available. The **ERA** handpiece can also use compatible universal cartridges.

18. GENERAL ADVICES FOR USE

COLOR TAKE

Dip the tip of the cartridge in a color holding capsule (tank), filled to the brim with color.

This operation will be carried out with the device turned on.

Through the Venturi effect the color tank fills up, allowing the needle to deposit pigment into the skin with each stroke.

ADVICE

Remember to periodically dip the tank into the color capsule, to avoid the needle working dry and therefore without adding pigment to the area to be treated.

19. TROUBLESHOOTING

THE DEVICE TURN OFF SUDDENLY

- Check the connection to the electrical network.
- Check that there is voltage in the network.
- Check that the cartridge is screwed in tightly.
- Check the integrity of the connection cable.
- Check the correct insertion of the connectors.

THE DEVICE HAS ABNORMAL VIBRATIONS

- Check the correct fitting of the cartridge.
- Check that the engine has not suffered impacts or decentralization.
- Check that the needle slides freely in the tube, without impediments.

THE DEVICE HEATS

- Check that the engine has not suffered impacts or decentralization.
- Check that the voltage is correct.
- Check that the accessories are mounted correctly.

THE DEVICE ONLY WORKS WITHOUT CARTRIDGE

- Try turning the device off and on again.
- Try replacing the cartridge.

IF IN DOUBT, IMMEDIATELY CONTACT THE BIOTEK ASSISTANCE CENTER NEAREST YOU OR THE BIOTEK ASSISTANCE LABORATORY IN MILAN DIRECTLY. IN ANY CASE SUSPEND ONGOING TREATMENTS.

20. MAINTENANCE AND CLEANING

The device must be maintained and used with all the precautions that a precision instrument requires.

Constant cleaning and disinfection of your device is of fundamental importance, to always offer your customer maximum safety.

For external cleaning it is necessary to use only specific products with soft, non-abrasive cloths.

Do not use excessively acidic or alkaline (aggressive) detergents which could compromise some components of the device.

To carry out the cleaning you need:

Disconnect the handpiece from the B4U Console or any power supply.

Disconnect external devices and other cables from the product.

Keep the product away from liquids.

Do not allow moisture to enter the openings and do not use sprays, solvents or abrasive substances.

Do not spray detergents directly on the product. Do not immerse in water or disinfectant solutions.

Do not use autoclaves or ultrasonic cleaners.

The cartridges must absolutely be replaced after each treatment and with each new customer.

Store the device in a closed, dust-free and aseptic place.

The device does not require any ordinary maintenance, other than cleaning and disinfection.

If the need for extraordinary interventions arises, contact only the BIOTEK laboratory in Milan.

For no reason is it permitted to intervene on the device arbitrarily or by personnel not authorized by the BIOTEK laboratory in Milan.

Carrying out a scheduled check-up every 12 months is absolutely recommended to keep your device efficient and safe.

21. DECLARATION OF CONFORMITY

BIOTEK s.r.l. declares the following:

The **ERA** device is CE marked and complies with the following directives:

- **ELECTROMAGNETIC COMPATIBILITY DIRECTIVE**
EMC Directive 2014/30/EC
- **LOW VOLTAGE DIRECTIVE**
LOW VOLTAGE Directive 2014/35/EC
- **MACHINERY DIRECTIVE**
MACHINERY Directive 2006/42/EC
- **FCC RULES FOR ELECTRONIC DEVICE - US**
CODE OF FEDERAL REGULATIONS, TITLE 47, PART 15 (47 CFR 15)

In compliance with the aforementioned directives, the following harmonized standards have been applied:

CEI EN 60335-1 Safety of domestic and similar devices
CEI EN 60335-1 Household and similar electrical devices safety

CEI EN 60601-1 Safety of electromedical device
CEI EN 60601-1 Medical electrical device safety – Part 1

CEI EN 55014-1 Electromagnetic compatibility
CEI EN 61000-3 Electromagnetic compatibility

In compliance with Interministerial Decree n° 206 of 15 October 2015 - Sheet n° 23

Certificate of origin

*This is a MADE IN ITALY product, designed and manufactured by BIOTEK with primary technology and quality materials. **ERA** by Biotek is the original Made in Italy.*

22. WARRANTY CERTIFICATE

This product was made with quality materials and careful specialized workmanship.
It was subjected to careful monitoring and testing by BIOTEK.

22. 1 BIOTEK device is guaranteed for 12 months from delivery directly by the manufacturer pursuant to art. 1495 Civil Code and according to the terms and conditions in force. The warranty gives the right to free repair of the product or parts of the product affected by manufacturing defects detected at the time of purchase or which appear during the warranty period.

Damages and breakages occurring during transport, damage resulting from tampering, repairs by third parties, improper use and maintenance are not covered by the warranty. You will also lose the warranty if liquids enter the engine because you have used non-original Biotek cartridges and accessories. We can guarantee a well-functioning workflow only if you use original BIOTEK products or compatible ones.

Materials subject to wear over time, such as the connection cable, connector and rubber O-rings, are also excluded from the warranty.

Compensation for any damage to people and/or things resulting from the use or failure of the product is excluded.

22. 2 **Under warranty.** In the unfortunate event that a Biotek machine, purchased on the Biotek shop or at one of the Biotek direct offices, had a malfunction/manufacturing defect manifested within the warranty period (12 months from the date of purchase), it must be shipped at care and expenses by the customer to BIOTEK in Milan for the repair. Under no circumstances will it be possible to issue a refund to the customer if the deadline for requesting a refund has already passed (see point 5). The costs of returning the device (unless a non-compliant use of the same is identified) will be borne by BIOTEK. Consult the page with the user manual for your machine to find out what your warranty includes.

22. 3 **Out of warranty.** Technical assistance interventions outside the warranty period require shipment to BIOTEK by the customer. The costs of return and technical interventions (including fault finding costs) will be charged to the customer and must be paid before processing.

22. 4 For any questions regarding device assistance, the customer can contact BIOTEK assistance on +39.02 7380144/145

23. RETURNS

23. 1 For purchases of device and/or accessories made via the internet, e-mail or telephone and in any case outside the commercial premises, it will always be possible to exercise the right of withdrawal within 14 calendar days from delivery of the goods (or the last good for multiple orders), contacting customer service for the return methods. The goods must be returned at the consumer's expense within 14 calendar days from the date on which the consumer communicated to BIOTEK his decision to withdraw.

23. 2 For purchases made through online or physical retailers other than Biotek srl, the customer is required to contact the company from which he purchased the goods directly. Biotek srl will not accept returns, replacements or repairs from customers who have not purchased directly from Biotek.

23. 3 IF the order was placed directly through the Biotek online shop and the return is requested within the technical times set out in point 5.1, BIOTEK will accept the return only and exclusively if the returned goods are in perfect condition and in their original packaging, accompanied by the relevant purchase invoice. In the event that the returns have been made in accordance with what is described above, BIOTEK will issue a refund via PayPal service or bank transfer according to the payment method made by the customer at the time of purchase, within a maximum of 14 days from receipt of the return at the BIOTEK headquarters in Milan. BIOTEK will not refund either the shipping costs invoiced with the purchase, nor the return costs.

23. 4 However, returns of products purchased directly from BIOTEK offices will not normally be accepted. In particular, once the BIOTEK sales office has been left, it will no longer be possible to exchange or replace the product: it is therefore up to the customer to verify, at the time of delivery, the quantity and condition of the products and, in the event of anomalies and/or of missing items/products, report it to BIOTEK staff before leaving the premises. An exception is made in cases in which the anomalies cannot be detected by the customer at the time of purchase: in this case replacements will be accepted, subject to verification by BIOTEK.

24. DISPUTES

The Court of Milan has exclusive jurisdiction for any dispute.

Biotek srl, via R. Farneti 8, 20129 Milano Italy
Tel. +39 (0)2.73.80.144/45 - Fax +39 (0)2.73.80.147 - info@biotek.it

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